

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



MCE MedChemExpress

Product Data Sheet

Derazantinib

 Cat. No.:
 HY-19981

 CAS No.:
 1234356-69-4

 Molecular Formula:
 C29H29FN4O

Molecular Weight: 468.57

Target: FGFR

Pathway: Protein Tyrosine Kinase/RTK

Storage: Powder -20°C 3 years 4°C 2 years

In solvent -80°C 2 years

-20°C 1 year

SOLVENT & SOLUBILITY

In Vitro

DMSO: 25 mg/mL (53.35 mM; Need ultrasonic)

Preparing Stock Solutions	Solvent Mass Concentration	1 mg	5 mg	10 mg
	1 mM	2.1342 mL	10.6708 mL	21.3415 mL
	5 mM	0.4268 mL	2.1342 mL	4.2683 mL
	10 mM	0.2134 mL	1.0671 mL	2.1342 mL

Please refer to the solubility information to select the appropriate solvent.

In Vivo

- 1. Add each solvent one by one: 10% DMSO >> 40% PEG300 >> 5% Tween-80 >> 45% saline Solubility: ≥ 2.5 mg/mL (5.34 mM); Clear solution
- Add each solvent one by one: 10% DMSO >> 90% corn oil Solubility: ≥ 2.5 mg/mL (5.34 mM); Clear solution

BIOLOGICAL ACTIVITY

Description

Derazantinib (ARQ-087) is an ATP competitive and orally active FGFR inhibitor (IC₅₀s: 1.8 nM for FGFR2, 4.5 nM for FGFR1 and 3 nM). Derazantinib inhibits FGFR phosphorylation. Derazantinib inhibits tumor growth in multiple xenograft models^{[1][2]}.

In Cells, inhibition of FGFR2 auto-phosphorylation and other proteins downstream in the FGFR pathway (FRS2α, AKT, ERK) is evident by the response to Derazantinib treatment. Cell proliferation studies demonstrate Derazantinib has anti-

proliferative activity in cell lines driven by FGFR dysregulation, including amplifications, fusions, and mutations. Cell cycle studies in cell lines with high levels of FGFR2 protein show a positive relationship between Derazantinib induced G1 cell

cycle arrest and subsequent induction of apoptosis^[1]. Derazantinib rescues the FGF2-mediated growth arrest with EC₅₀ at about 100 nM, with no significant toxicity detected for up to 500 nM. The concentration range at which Derazantinib significantly suppresses the FGF2 effect is between 70-500 nM. Derazantinib inhibits FGF-mediated loss of extracellular matrix and induction of chondrocyte premature senescence. Derazantinib rescues FGF-mediated inhibition of chondrocyte differentiation in tibia cultures. Derazantinib inhibits FGFR1-4 but no other receptor tyrosine kinases in cell-free kinase assay. Derazantinib inhibits FGFR1 and FGFR2 mutants associated with craniosynostoses. Derazantinib rescues FGFR-mediated bone differentiation in mouse limb bud micromass cultures and ex vivo mouse calvarial organ cultures^[2]. MCE has not independently confirmed the accuracy of these methods. They are for reference only.

In Vivo

Derazantinib is effective at inhibiting tumor growth in FGFR2 altered, SNU-16 and NCI-H716, xenograft tumor models with gene amplifications and fusions^[1]. Most of the embryos exhibit abnormal external phenotype (81.3%) in Derazantinibinjected wings, possibly due to inhibition of proliferation of limb bud mesenchyme. The wings are shorter and thinner, with skeletal phenotype typical for FGFR inhibition, where ulna and radius are shorter or smaller in size, or occasionally missing completely^[2].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

PROTOCOL

Kinase Assay [1]

Derazantinib is titrated in DMSO utilizing a 3-fold dilution scheme, and then diluted 10-fold further in deionized water for a final DMSO concentration of 10%. A volume (2.5 μ L) of these dilutions or vehicle is added to each well of a reaction plate. FGFR1 or FGFR2 is added to assay buffer to each well in a volume of 17.5 μ L for a final concentration of 0.50 or 0.25 nM, respectively. After a 30-minute pre-incubation period, ATP and substrate are added in assay buffer (5 μ L) for final concentrations of 0-1,000 μ M ATP and 80 nM biotinylated-PYK2, for a final reaction volume of 25 μ L. The plates are incubated for 60 minutes at room temperature, and then stopped in the dark by the addition of 10 μ L stop/detection mixture prepared in assay buffer containing EDTA^[1].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

Cell Assay [1]

Cells are seeded at 3000-5000 cells per well with 130 μ L media in 96-well tissue culture treated plates. The cells are incubated overnight and subsequently treated with 3-fold serial dilutions of Derazantinib starting at 100 μ M. The cells are returned to a 37°C humidified incubator for 72 hours. Cell proliferation is measured using MTS assay^[1].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

Animal
Administration [1]

Mice: Female NCr nu/nu mice (SNU-16) or CB17 SCID mice (NCI-H716) with well established (400 mg) subcutaneous tumors are given a single oral dose of Derazantinib or vehicle control. Plasma and tumor samples are collected 4 hours post single dose. Derazantinib is administered orally. Dosing volume for all groups is 10 mL/kg or 0.1 mL/10 g body weight [1].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

CUSTOMER VALIDATION

• Patent. US11768205B2.

See more customer validations on www.MedChemExpress.com

REFERENCES

[1]. Hall TG, et al. Preclinical Activity of ARQ-087, a Novel Inhibitor Targeting FGFR Dysregulation. PLoS One. 2016 Sep 14;11(9):e0162594.

[2]. Balek L, et al. ARQ-087 inhibits FGFR signaling and rescues aberrant cell proliferation and differentiation in experimental models of craniosynostoses and

		F0FD0 F0FD0 D 00175	105.57.00	
nondrodysplasias caused by a	ctivating mutations in FGFR1,	FGFR2 and FGFR3. Bone. 2017 [Jec;105:57-66.	
	Cautian, Draduct has not	hoon fully validated for mo	dieal anniiestiana Fayyaaayah	مام
			dical applications. For research use o	
	Tel: 609-228-6898	Fax: 609-228-5909	E-mail: tech@MedChemExpress.	com
	Address: 1 D	eer Park Dr, Suite Q, Monmo	uth Junction, NJ 08852, USA	

Page 3 of 3 www.MedChemExpress.com