

Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com



www.rockland.com tech@rockland.com +1 484.791.3823

Datasheet for R109-0050 Rabbit Blood Defibrinated

Overview

Description:	Rabbit Blood Defibrinated - R109-0050
Item No.:	R109-0050
Size:	50 mL
Applications:	Other
Origin:	Rabbit

Product Details

Background:	Defibrinated rabbit blood is aseptically collected and prepared to order from our colony of animals. Fibrin is mechanically removed during the clotting process without the presence of anticoagulants or other additives. Defibrinated rabbit blood is ideal for blood agar plate preparations. Blood Agar is a general purpose enriched medium often used to grow fastidious organisms and to differentiate bacteria based on their hemolytic properties. Defibrinated rabbit blood is perishable and is collected and processed upon receipt of your order.
Synonyms:	Sterile Rabbit Blood Defibrinated
Species of Origin:	Rabbit

Target Details

Relevant Links:	•	R109 SDS

Application Details

Suggested Applications:	Other (Based on references)
Application Note:	Blood agar plate preparations
Assay Dilutions:	All assays should be optimized by the user. Recommended dilutions (if any) may be listed below.

Tissue Data



www.rockland.com tech@rockland.com +1 484.791.3823

Tissue Type:	Blood
Sex:	Mixed
Strain:	Rabbit - New Zealand White


Formulation

Physical State:	Liquid
Buffer:	None
Anticoagulant:	Defibrinated
Sterility:	Sterile
Preservative:	None
Stabilizer:	None

Shipping & Handling

Shipping Condition:	Wet Ice
Storage Condition:	Store defibrinated rabbit blood at 4° C prior to opening. Use aseptic technique to maintain sterility when opening product. Be advised that blood is a perishable product and exact shelf may depend on application.
Expiration:	This product MAY be stable for up to two (2) weeks if properly stored and handled.

Images



ELISA

Comparison of blood cyanide concentrations measured at UCSD and UTA by two different methods. Experimental rabbit at UCI was infused with 0.33 mg/min cyanide and samples taken at different times. Defibrinated rabbit blood (p/n R109-0100). Figure 4. PMID: 23473259.



www.rockland.com tech@rockland.com +1 484.791.3823

References

- Singh N, Wang C, Cooper R. Potential of Essential Oil-Based Pesticides and Detergents for Bed Bug Control. *J Econ Entomol.* (2014)
- Tian Y et al. A disposable blood cyanide sensor. Anal Chim Acta. (2013)

Disclaimer

This product is for research use only and is not intended for therapeutic or diagnostic applications. Please contact a technical service representative for more information. All products of animal origin manufactured by Rockland Immunochemicals are derived from starting materials of North American origin. Collection was performed in United States Department of Agriculture (USDA) inspected facilities and all materials have been inspected and certified to be free of disease and suitable for exportation. All properties listed are typical characteristics and are not specifications. All suggestions and data are offered in good faith but without guarantee as conditions and methods of use of our products are beyond our control. All claims must be made within 30 days following the date of delivery. The prospective user must determine the suitability of our materials before adopting them on a commercial scale. Suggested uses of our products are not recommendations to use our products in violation of any patent or as a license under any patent of Rockland Immunochemicals, Inc. If you require a commercial license to use this material and do not have one, then return this material, unopened to: Rockland Inc., P.O. BOX 5199, Limerick, Pennsylvania, USA.