

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



PI 4-kinase α siRNA (h): sc-44012



The Power to Question

BACKGROUND

The members of the phosphatidylinositol kinase (PIK) superfamily can be divided into three groups based on their substrate specificity. PIKs convert phosphatidylinositol (PI) into PI phosphate [PI(3)P], PI phosphate [PI(4)P], PI bisphosphate [PI(4,5)P2] and PI triphosphate [PI(3,4,5)P3]. The first group, the PI 3-kinases, is composed of highly related proteins designated p110 α , p110 β , p110 γ and p110 δ which convert PI into PI(3)P and PI(4,5)P2 into PI(3,4,5)P3. The second group, the PI 4-kinases, convert PI into PI(4)P. The third group, the PI(4)P5-kinases, convert PI(4)P into PI(4,5)P2. Phosphatidylinositides represent important regulatory molecules and are involved in a diverse array of signaling pathways. Phosphatidylinositol biphosphate acts as an activator of PKCs and as a substrate for PLC γ , which converts the molecule into the second messengers, inositol-1,4,5 triphosphate and 1,2-diacylglycerol. PI(3,4,5)P3 has been shown to activate the PKC ζ isoform. Wortmannin, originally described as a specific inhibitor of PI 3-kinases, may actually be a broad spectrum inhibitor of PI kinase activity.

REFERENCES

- Hara, K., et al. 1994. 1-phosphatidylinositol 3-kinase activity is required for Insulin-stimulated glucose transport but not for RAS activation in CHO cells. Proc. Natl. Acad. Sci. USA 91: 7415-7419.
- Roche, S., et al. 1994. The phosphatidylinositol 3-kinase a is required for DNA synthesis induced by some, but not all, growth factors. Proc. Natl. Acad. Sci. USA 91: 9185-9189.
- Stephens, L., et al. 1994. A novel phosphatidylinositol 3 kinase activity in myeloid-derived cells is activated by G protein βγ subunits. Cell 77: 83-93.
- 4. Woscholski, R., et al. 1994. Biochemical characterization of the free catalytic p110 α and the complexed heterodimeric p110 α .p85 α forms of the mammalian phosphatidylinositol 3-kinase. J. Biol. Chem. 269: 25067-25072.
- Woscholski, R., et al. 1994. A comparison of demethoxyviridin and wortmannin as inhibitors of phosphatidylinositol 3-kinase. FEBS Lett. 342: 109-114.
- 6. Hunter, T. 1995. When is a lipid kinase not a lipid kinase? When it is a protein kinase. Cell 83: 1-4.

CHROMOSOMAL LOCATION

Genetic locus: PI4KA (human) mapping to 22q11.21.

PRODUCT

PI 4-kinase α siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see PI 4-kinase α shRNA Plasmid (h): sc-44012-SH and PI 4-kinase α shRNA (h) Lentiviral Particles: sc-44012-V as alternate gene silencing products.

For independent verification of PI 4-kinase α (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-44012A, sc-44012B and sc-44012C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNAse-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNAse-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

PI 4-kinase α siRNA (h) is recommended for the inhibition of PI 4-kinase α expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 µM in 66 µl. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor PI 4-kinase α gene expression knockdown using RT-PCR Primer: PI 4-kinase α (h)-PR: sc-44012-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

Santa Cruz Biotechnology, Inc. 1.800.457.3801 831.457.3801 fax 831.457.3801 Europe +00800 4573 8000 49 6221 4503 0 www.scbt.com