



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

# Myozenin 1 siRNA (m): sc-149767

## BACKGROUND

Myozenins, including Myozenin 1 (also known as MYOZ1, Calsarcin-2, or FATZ), Myozenin 2 and Myozenin 3, are a family of intracellular binding proteins that are involved in linking Z-disk proteins to the sarcomere. Expressed at high levels in cardiac and skeletal muscle and at lower levels in prostate, pancreas and heart, Myozenin 1 is a 299 amino acid nuclear protein that interacts with the Z-disk proteins,  $\alpha$ -actinin, Filamin 2 and PP2B (calcineurin), and effectively forms a bridge between these proteins and muscle fibers. Via these interactions, Myozenin 1 couples striated muscle activity to protein activation and is thought to play a role in both Z-disk assembly and myofibrillogenesis. Due to the close involvement of Myozenin 1 with muscle formation, mutations in the gene encoding Myozenin 1 may be associated with muscular dystrophies and neuromuscular myopathies.

## REFERENCES

1. Faulkner, G., et al. 2000. FATZ, a filamin-, actinin-, and telethonin-binding protein of the Z-disc of skeletal muscle. *J. Biol. Chem.* 275: 41234-41242.
2. Frey, N., et al. 2000. Calsarcins, a novel family of sarcomeric calcineurin-binding proteins. *Proc. Natl. Acad. Sci. USA* 97: 14632-14637.
3. Takada, F., et al. 2001. Myozenin: an  $\alpha$ -actinin- and  $\gamma$ -filamin-binding protein of skeletal muscle Z lines. *Proc. Natl. Acad. Sci. USA* 98: 1595-1600.
4. Online Mendelian Inheritance in Man, OMIM<sup>™</sup>. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 605603. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
5. Gontier, Y., et al. 2005. The Z-disc proteins myotilin and FATZ-1 interact with each other and are connected to the sarcolemma via muscle-specific filamins. *J. Cell Sci.* 118: 3739-3749.
6. Arola, A.M., et al. 2007. Mutations in PDLIM3 and MYOZ1 encoding myocyte Z line proteins are infrequently found in idiopathic dilated cardiomyopathy. *Mol. Genet. Metab.* 90: 435-440.
7. Posch, M.G., et al. 2007. Mutations in MYOZ1 as well as MYOZ2 encoding the calsarcins are not associated with idiopathic and familial dilated cardiomyopathy. *Mol. Genet. Metab.* 91: 207-208.

## CHROMOSOMAL LOCATION

Genetic locus: Myoz1 (mouse) mapping to 14 A3.

## PRODUCT

Myozenin 1 siRNA (m) is a pool of 2 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Myozenin 1 shRNA Plasmid (m): sc-149767-SH and Myozenin 1 shRNA (m) Lentiviral Particles: sc-149767-V as alternate gene silencing products.

For independent verification of Myozenin 1 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-149767A and sc-149767B.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

Myozenin 1 siRNA (m) is recommended for the inhibition of Myozenin 1 expression in mouse cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## GENE EXPRESSION MONITORING

Myozenin 1 (8): sc-136429 is recommended as a control antibody for monitoring of Myozenin 1 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG $\kappa$  BP-HRP: sc-516102 or m-IgG $\kappa$  BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker<sup>™</sup> Molecular Weight Standards: sc-2035, UltraCruz<sup>®</sup> Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG $\kappa$  BP-FITC: sc-516140 or m-IgG $\kappa$  BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz<sup>®</sup> Mounting Medium: sc-24941 or UltraCruz<sup>®</sup> Hard-set Mounting Medium: sc-359850.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Myozenin 1 gene expression knockdown using RT-PCR Primer: Myozenin 1 (m)-PR: sc-149767-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.