



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

P.O. Box 3286 - Logan, Utah 84323, U.S.A. - Tel. (800) 729-8350 - Fax (435) 755-0015 - www.scytek.com

HRP Stabilizing Diluent

Description: HRP Stabilizing Diluent (MOPS Buffered) has been developed specifically for the stabilization of Horseradish Peroxidase conjugates. HRP Stabilizer prevents the loss of catalytic activity and maintains the conformation of the antibody or protein antigen portion of the conjugate. Product is filtered at 0.2 microns.

Contents: Proprietary mixture of HRP stabilizers and enhancers.

Preservatives: 0.02% Bromonitrodioxane, 0.02% 2-Methyl-4-isothiazolin-3-one, 0.005% Proclin 300 (Proclin is a registered trademark of Rohm & Haas Company), and 0.0008% Neomycin.

Uses/Limitations: For dilution and stabilization of HRP conjugates in solution for immunochemical procedures.
Do not use past expiration date.
Do not use if solution becomes cloudy.
Blending with phosphate buffered saline may result in precipitation. To prevent precipitation, use normal saline (0.85% NaCl) rather than phosphates.
Product does not alter the stability of the bond between the enzyme and the antibody or antigen.

Availability:


<u>Item #</u>	<u>Volume</u>
HSM125	125 ml
HSM500	500 ml
HSM999	1000 ml


Please contact for bulk and OEM volumes and pricing.

Storage: Store at 2-8°C. Solution is stable for at least 24 months from date of manufacture. Significantly increases shelf-life of HRP once diluted into this solution. Product may be stored room temperature for short periods (1 – 2 weeks). May be shipped at ambient temperature or with wet ice. Avoid continuous exposure to direct light. Do not freeze.


Procedure: Dilute HRP to concentration required for the assay. Store diluted solution at 2-8°

Precautions: Avoid contact with skin and eyes.
Wash after use.
Observe all federal, state and local environmental regulations regarding disposal.

Storage: 2° C  8° C

 ScyTek Laboratories, Inc.
205 South 600 West
Logan, UT 84321
435-755-9848
U.S.A.


Emergo Europe
Prinsessegracht 20
2514 AP The Hague, The Netherlands