



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Azlocillin sodium salt

Chemical Properties


CAS No. : 37091-65-9

Formula: C₂₀H₂₂N₅NaO₆S

Molecular Weight: 483.47

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Azlocillin sodium, an acylampicillin, is a broad-spectrum antibiotics
Targets(IC50)	Antibacterial,Antibiotic,Parasite
In vitro	Azlocillin (12.5 µg/mL) inhibits over 75% of the isolates of <i>Pseudomonas aeruginosa</i>. Azlocillin (12.5 µg/mL) is also active against indole-negative and -positive <i>Proteus</i> spp., inhibiting 98 and 71%, respectively. Azlocillin is more active than mezlocillin, ticarcillin, and carbenicillin and as active as BLP-1654 against isolates of <i>P. aeruginosa</i>. [1] The acyl side chains of Azlocillin have an ureido-(urea) structure hence the name 'ureidopenicillins' or, more specifically, 'acylureidopenicillins.' In vitro studies against <i>P. aeruginosa</i> demonstrates that piperacillin has activity that is twice that of azlocillin, 4 times that of mezlocillin and ticarcillin, and about 8 times that of carbenicillin. Azlocillin produces elongated bacterial forms with delayed or no lysis in morphologic studies. [2] Azlocillin has MICs of 12.5 µg/mL on <i>Pseudomonas aeruginosa</i>. Azlocillin (3.125 µg/mL) results in a reduction in the rate of growth but no bactericidal phase on <i>Pseudomonas aeruginosa</i>. Azlocillin decreases an initial lag phase with increasing drug concentration. At the lower concentration of tobramycin (0.5 µg/ml), the combinations with both the high and the low concentrations of Azlocillin are more effective than the individual components on <i>Pseudomonas aeruginosa</i>. [3] Isolates with derepression of AmpC enzyme are one to two doubling dilutions more resistant to azlocillin than are those in which increased efflux or impermeability is inferred. Those with secondary β-lactamases are mostly (12/14 cases) susceptible to ceftazidime at 4 mg/L, but are amongst the most resistant to Azlocillin (MIC ≥128 mg/L in 10/14 cases). [4]

Solubility Information

Solubility	DMSO: 90 mg/mL (186.2 mM), H ₂ O: 89 mg/mL (184.1 mM), Ethanol: 3 mg/mL (6.19 mM), (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

A DRUG SCREENING EXPERT

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0684 mL	10.3419 mL	20.6838 mL
5 mM	0.4137 mL	2.0684 mL	4.1368 mL
10 mM	0.2068 mL	1.0342 mL	2.0684 mL
50 mM	0.0414 mL	0.2068 mL	0.4137 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Stewart D, et al. Antimicrob Agents Chemother, 1977, 11(5), 865-870.

Wright AJ, et al. Mayo Clin Proc, 1999, 74(3), 290-307.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481